

ReSound GN

Next-era hearing

Top-rated for hearing in noise.
Smaller than ever.
Uncompromised. And...

ReSound Nexia



Smaller than ever. Uncompromised. And...

Introducing ReSound Nexia™, a ground-breaking sound experience for mild to profound hearing loss. Top-rated for hearing in noise. Ready for Bluetooth® Auracast™. And... with our most advanced technology and features to support your clients.

... powered by Organic Hearing

Our unique approach to creating hearing solutions that sound natural, feel natural and connect your clients naturally to the world. Discover the difference now.



... connectivity with a human touch

Ready for Bluetooth Auracast, the future of public connectivity to be rolled out globally across all devices, giving users the power to connect to anything everywhere

- New TV-Streamer+ based on Bluetooth Auracast for clear audio at the preferred volume: Press once to pair with ReSound Nexia
- Swipe to access the new sound program and connect to the TV-Streamer+ in the ReSound Smart 3D™ app
- Stream calls with a double tap on the ear or the hearing aids



...top-rated for hearing in noise



89%

... recommended by users

Don't just take our word for it. 89% of users recommend ReSound Nexia technology for hearing in noise.¹ Discover the difference ReSound Nexia makes for your clients — both for newcomers to hearing aids or existing clients considering an upgrade.



+50%
user satisfaction²

+48%
fit acceptance³

... accepted by users

ReSound Nexia helps solve for the greatest pain points for hearing aid users — hearing in noise, unnoticeability, wearing comfort and ease-of-use. Because of this, we're seeing higher acceptance from our users than ever before, and higher overall satisfaction².



... a powerful miniaturised design without compromises

The rechargeable microRIE is our smallest style ever. Yet it's capable of an astonishing 150% improvement in speech understanding in noise³.

- Expertly crafted to ensure the best possible ergonomics and mechanics
- Sleek organic shape is 25% smaller than the standard RIE⁴
- All-day power from a short recharge
- IP68-rated all-weatherproof designs

³ 4.36 dB SNR improvement compared to legacy products. Jespersen & Groth (2022)

⁴ Technical data on file

... flexible fitting options for all needs

Give your clients a great first experience with hearing aids and enable them to hear their best in noise with ReSound Smart Fit™ fitting software.

- The new experience profile 'First Time User Onboarding' is ideal for initial sound comfort and adaptation to hearing amplification
- Give your clients the best fit for hearing in noise right from the start with the new 'Hear in Noise' program
- Offer a superb wireless CROS/BiCROS solution with greater audibility all around, with the new solution for single-sided deafness
- Perfectly complementing Cochlear™ implants



¹ Compared to legacy technology. End-user NPS rating, data on file
² Compared to legacy technology. End-user NPS rating, data on file
³ Compared to legacy technology. GN RFC data on file

Your style options



Key feature overview

Features	NX9	NX7	NX5	NX4
Sound Quality				
WARP bands	★★★★ 17	★★★★ 14	★★★☆☆ 12	★★★☆☆ 12
Environment recognition	★★★★ (Environmental Optimiser II)	★★★☆☆ (Environmental Optimiser)	★★★☆☆ (Environmental Optimiser)	
Directionality	★★★★ (360 All-Around)	★★★☆☆ (All-Access)	★★★☆☆ (Binaural Directionality III)	★★★☆☆ (Binaural Directionality)
Speech intelligibility in noise	★★★★ (Front Focus)	★★★☆☆ (Ultra Focus)	★★★☆☆ (Adaptive Directionality)	★★★☆☆ (Adaptive Directionality)
Noise Reduction	★★★★ (5 settings)	★★★☆☆ (3 settings)	★★★☆☆ (2 settings)	★★★☆☆ (On / Off)
Comfort	★★★★ (3 settings)	★★★★ (3 settings)	★★★☆☆ (On / Off)	★★★☆☆ (On / Off)
Feedback Management	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)
Low frequency boost	★★★★ (3 settings)	★★★☆☆ (2 settings)	★★★☆☆ (On / Off)	★★★☆☆ (On / Off)
Connectivity				
Direct audio streaming	★★★★★ For iPhone, iPad, iPod touch and Android™ devices. Hands-free calls for iPhone and iPad*			
Public connectivity	★★★★★ Future ready with Auracast for new broadcast opportunities, telecoil built into RIE 13 battery style and in Multi Mic			
Accessories	★★★★★ TV-Streamer+, TV Streamer 2, PhoneClip+, Multi Mic, Micro Mic, Remote Control, ReSound Smart 3D App			
Fitting Flexibility				
Wireless CROS	★★★★★			
Gain handles	17	14	12	12
Programs	4	4	4	4
Fitting options	Yes	Yes	Yes	Yes
Bimodal compatibility	★★★★★			

* Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.

Worldwide Headquarters

GN ReSound A/S
DK-2750 Ballerup
Denmark
Tel.: +45 4575 1111
resound.com

CVR no. 55082715

United Kingdom

GN Hearing UK Ltd.
Unit 13
Talisman Business Centre
Bicester
Oxon OX26 6HR
Tel: +44 1869 352800
resound.com

Australia

GN Hearing Australia Pty Ltd
Gate C, 19-25 Khartoum Road
Macquarie Technology Park
Macquarie Park NSW 2113
Tel.: (free) 1800 658 955
resound.com

New Zealand

GN Hearing New Zealand
Limited
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
resound.com